
Schmelzpunktbestimmungsapparate MPM-H3 und MPM-HV3
Melting point determination units MPM-H3 and MPM-HV3

OPERATING INSTRUCTIONS

(DEUTSCH / ENGLISH / FRANÇAIS / ESPAÑOL / РУССКИЙ)

Messung nach Arbeitsanweisung PH.- Eur. Methode
Measurement according to work instruction PH.- Eur. Method
Mesure selon les consignes de travail PH- Eur. Méthode
Medición de acuerdo con lo establecido en las instrucciones de trabajo Método de la Farmacopea Europea
Измерения согласно методам рабочей инструкции ЕФ

INDEX

Deutsch

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PRODUCT DATA

GENERAL INFORMATION

1. Only pull the safety hood over the device after cooling off.
2. Always insert the capillaries with the closed side first.

TECHNICAL DATA

Mains voltage wide-range input	88 – 264 V / AC 0.8 A
Heating output	12 V / DC 25 VA
Protection class	IP20
Measuring range	up to 360 °C
Measuring accuracy	+/- 0.4 °C
Reproducibility	+/- 0.3 °C
Heating rate	35 °C to 100 °C approx. 2 minutes 35 °C to 300 °C approx. 7 minutes
Cooling off	100 °C to 35 °C approx. 4 minutes 300 °C to 35 °C approx. 9 minutes
Magnifying glass	20 dioptres / 5-fold
Capillaries (outside x inside x length)	1.45 x 1.0 x 80 mm unilaterally closed.
Dimensions housing (W x H x D)	32 cm x 26 cm x 6.5 cm
Dimensions display (W x H)	8.5 cm x 5.5 cm
Weight of device	4 kg
Weight of stand	2.8 kg
Height adjustment	From table top to middle of magnifying glass from 36 cm to 46 cm

ACCESSORIES

- 1x mains cable
- 1x safety hood
- 1x 100 capillaries
- 1x USB cable
- 1x software as CD or USB stick
- 1x operating instructions

SAFETY NOTES / COMMISSIONING

ATTENTION! SAFETY NOTES

- Do not cover the ventilation slits
- Do not push any objects into the ventilation slits
- No moisture must enter the device
- Only pull the safety hood over the device after cooling off.
- The capillaries will heat more at several subsequent melting point measurements
- Only use the intended capillaries for melting point determinations
- Observe the corresponding safety provisions when handling chemicals
- Always disconnect the unit from the mains by pulling the mains plug before any maintenance and care work at the device.
- Work and servicing at the device must only be performed by authorized specialists or the manufacturer
- Take the unit off of the stand for servicing and pack it safely for transport without accessories

UNPACK – INSTALLATION

- Unpack the unit with care
- Report transport damage at once
- Check completeness according to the accessories list
- Let the device adjust to room temperature

COMMISSIONING

Find a suitable location for the device. Avoid extreme proximity to heat sources.

LANGUAGE SETTINGS

Keep the set-button pushed and switch on the device. Use an arrow-button (temperature preselection) to select the operating language and confirm it with the start-button.

SETTING THE TIME AND DATE

1. First switch on the device, then push the menu- and set-buttons at the same time.
2. Use the arrow-buttons to set the date and day. Then use the start-button to switch to the month and make your selection again with the arrow-buttons.
3. Follow these steps further until you have set the date and time.
4. Finally, leave the menu by pushing the start-button.

INSTRUCTIONS FOR THE PC SOFTWARE

If you have not installed any PDF converter on your computer, download a freely available version from the internet so that you can choose the file format PDF as a virtual printer.

INSTALLATION

AUTOMATIC INSTALLATION


Insert the CD in your computer. Start installer.exe and follow the installation instructions.

MANUAL INSTALLATION

Driver Installation under Windows XP

You may need administration rights in order to install the driver.

1. Switch on your MPM unit and connect the USB cable. The following window appears:

2. Click that you do not want to connect to Windows update and click **Next**


3. Make sure that the driver CD is in the CD/DVD drive. Select that you do not want to install the driver automatically and click **Next**

4. Select the Driver folder on the driver CD and click **Next**

5. The driver is installed. After successful installation, the following window appears:

6. Confirm with **Finish**.
You can use your MPM unit now with the MPM software.

MODEL MPM-H3


Determination of a melting point (visual)

1. Set the required starting temperature 3 °C below the expected melting point of your substance. Additionally pushing the set button during the temperature selection reduces the setting in steps of 10 °C.
2. Push the start button. The device starts to heat up to the set starting temperature. Once the temperature is reached, a sound signal will be emitted. The starting temperature in the display flashes. If no measurement is started, the device will automatically cool off after 12 minutes.
3. Insert the filled capillary (approx. 4 mm) into the intended capillary guide with the closed side first at the top of the device.
4. Confirm the start button (the magnifying glass lighting switches on automatically and the device switches to heating mode 1°C/Min). Observe the sample through the magnifying glass. When the melting point is reached, push the stop button. The melting point is displayed. If no melting point has been found after 12 minutes, the device will automatically switch to cooling, and a new measurement can be performed.
5. You can perform three measurements in sequence before the display shows

If you do not need the measured data, push the start button to delete the data. For a printout of the logs, start the program MPM on your computer.

Clicking the Download symbol will transfer the measuring results via the melting point determination unit to the computer.

The USB symbol must appear when the USB connection to the meter is active


6. There are three menu tabs above the figure with time indications of the times for the respective measurements. Clicking can query the individual measuring results.
7. Enter the necessary data for the respective measuring result in the mandatory fields.
 - Name of the person responsible for the measurement
 - Place where the measurement has been performed, e.g. lab
 - Date of the measurement
 - Name of the tested substance
 - Melting range from / to in °C of the substance according to the certificate
 - Batch or lot number of the substance
 - Expiration date of the substance
 - Field for additional comments
8. Use the printer icon/printer selection to choose SCAN soft PDF create. Then you can print out your measuring log as a PDF file or save it on your computer.

MODEL MPM-HV3

1. Set the required starting temperature 3 °C below the expected melting point of your substance. Additionally pushing the set button during the temperature selection reduces the setting in steps of 10 °C.
2. Push the start button. The unit starts to heat up to the chosen starting temperature. Once the temperature is reached, a signal sounds. The starting temperature in the display flashes. If no measurement is started, the device will automatically cool off after 12 minutes.
3. Choose the measuring procedure visual or automatic. In case of visual measurement through the magnifying glass continue as described in items 1 to 8, model MPM-H3.
4. For automatic measurement, insert the filled capillaries (approx. 4 mm) into the rear marked measuring channel (autom. \ominus) with the closed side first (First tap the capillaries down to compress the substance).
5. The measurement will start automatically at 1 °C / minute when inserting the capillaries.
6. During the measurement, you can track the change to the substance graphically illustrated in % on the display. In rare cases, bubbles may form in the capillary during melting, which may prevent detection of the melting point. In that case, the measurement must be repeated. If the device does not recognize any melting point after 12 minutes, the fan cooling will be activated automatically.
7. When the melting point is reached, a signal will sound and the display will indicate this. You can now download the data from a computer and save them as a PDF file and print them. After three subsequent measurements, the display will show

If you do not need the measured data, push the start button to delete the data. For a printout of the logs, start the program MPM on your computer.

Clicking the Download symbol will transfer the measuring results via the melting point determination unit to the computer.

The USB symbol must appear when the USB connection to the meter is active

8. There are three menu tabs above the figure with time indications of the times for the respective measurements. Clicking can query the individual measuring results.
9. Enter the necessary data for the respective measuring result in the mandatory fields.
 - Name of the person responsible for the measurement
 - Place where the measurement has been performed, e.g. lab
 - Date of the measurement
 - Name of the tested substance
 - Melting range from / to in °C of the substance according to the certificate
 - Batch or lot number of the substance
 - Expiration date of the substance
 - Field for additional comments
10. Use the printer icon/printer selection to choose SCAN soft PDF create. Then you can print out your measuring log as a PDF file or save it on your computer.

INSPECTION LOG

Example log with figure for an automatic measurement (MPM-HV3)


No figure will be created when measuring in the semi-automatic mode (visual).

After the log has been printed out, the control window below the respective log for printout below the figure will be checked automatically to indicate which logs have already been printed. If you want to print a measuring result again, click to deactivate the „printed“ symbol.

The measuring log will automatically show the serial number of the device so that the log can be assigned to the device according to QMS management.

To check the QMS management, the next due calibration date is indicated automatically at the end of the log.

The measuring log can be printed in the respective language set by the user, German or English.

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