---

---

# Dispenser / Dispenser / Distributeur / Dosificadore

Gebrauchsanweisung  
Operating Manual  
Mode d'emploi  
Instrucciones de manejo

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Before using the instrument for the first time, ensure it is rinsed carefully or discard the first few samples dispensed.

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## Contents

1. Safety Instructions	4
2. Application and Operating Limitations	5
3. Operating Exclusions	5
4. Operating Limitations	6
5. Components	7
6. Dispensing	8
6.1. Preparation for dispensing	8
6.2. Priming	9
6.3. Setting the volume	9
6.4. Dispensing	9
7. Cleaning	10
7.1. Emptying	10
7.2. Standard cleaning	11
7.3. Intensive cleaning	11
8. Cleaning/replacing valves	12
8.1. Filling valve	12
8.2. Discharge valve	12
9. Autoclaving	13
10. Volume Check	13
11. Troubleshooting	14
12. Repair Service	14
13. Warranty	14
14. Technical Data/Ordering Information	15
14.1. Items supplied	15
14.2. Accessories and Spare Parts	16

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# 1. Safety Instructions

This Manual does not purport to address every safety issue which may arise during use. It is the responsibility of the end user of this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Use the utmost caution when dispensing caustic, poisonous, radioactive or hazardous chemicals.
2. Observe general safety regulations (e.g., wear protective clothing, goggles and gloves).
3. Observe the Operating Manual and information from reagent manufacturers.
4. Use the instrument only for dispensing liquids, with strict regard to the defined Operating Exclusions and Limitations. If in doubt regarding the suitability of the instrument, for a particular application contact the manufacturer.
5. Before use, always verify that the instrument is in good working order, e.g., piston moves smoothly; filling and discharge tube and secure discharge valve are firmly seated and are properly attached, etc.
6. Never use force on the instrument. Use of force may result in injury to the user and other persons.
7. While dispensing, the discharge tube must always point away from the user or other persons. Avoid splashes. Only dispense into suitable vessels.
8. Never carry the mounted instrument by its upper casing.
9. Only operate the instrument with correctly installed discharge tube. Do not push down the piston while the closure cap is screwed on.
10. Clean the instrument before disassembling.
11. Use only original manufacturer's accessories and spare parts. Do not carry out any technical alterations.
12. In case of difficulty (e.g., piston difficult to move, leakage), immediately stop dispensing. Before further use, repair the instrument as described in this Manual. Contact the manufacturer if necessary.

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## 2. Application and Operating Limitations

This instrument is designed for dispensing liquids, observing the following physical limits:

- 15 to 40 °C for instrument and reagent
- density up to 2.2 g/cm<sup>3</sup>
- vapor pressure up 500 mbar
- viscosity up to 500 mm<sup>2</sup>/s

## 3. Operating Exclusions

When the instrument is correctly used, the dispensed liquid comes into contact with only the following materials:

Borosilicate glass, FEP, ETFE, PFA and Platinum-Iridium.

The instrument is suitable for dispensing a wide range of liquids, **excluding**:

- Liquids which attack FEP, ETFE and PFA;
- Solutions containing hydrofluoric acid, since they attack borosilicate glass;
- Suspensions, containing solid particles;
- Solutions which decompose and form solid particles (e.g., Biuret reagent);
- Substances which undergo catalytic transformation or react with Platinum-Iridium (e.g., H<sub>2</sub>O<sub>2</sub>);
- Carbon disulfide, as this media inflames easily

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## 4. Operating Limitations

Liquids which form deposits may make the piston difficult to move or may cause jamming (e.g., crystallizing solutions or highly concentrated alkaline solutions).

When dispensing flammable media, make provisions to avoid static charging, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.

### Warning!


If there is a sign of a potential malfunction (e.g., piston difficult to move) never use force. Immediately stop dispensing and follow cleaning instructions (see chapter 7) or contact the manufacturer.

### Note:

Compatibility of the instrument for this special application (e.g., trace material analysis) must be checked by the user or contact the manufacturer.

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## 6. Dispensing

### 6.1. Preparation for dispensing

1. Observe Operating Limitations and general safety regulations.

#### **Warning!**

Only the slotted-type discharge tube supplied with the instrument allows correct operation. Do not use any other type of discharge tube. Never use damaged or deformed discharge tubes.

2. Push the discharge tube as far as possible into the discharge valve.
3. Firmly finger-tighten the locking nut. Check the discharge tube for a tight fit.
4. The opening discharge tube must always point downwards.
5. Cut the length of the intake tube (according to the bottle to be used) and push the intake tube into the intake valve as far as it will go. Cut the lower end at an angle.
6. Mount the dispenser upon the reagent bottle.  
For bottles with different thread diameter select a suitable adapter (see chapter 14.2).
7. Hold the discharge tube support and carefully remove the discharge tube closure cap.

#### **Warning!**

Avoid splashing of reagent!

8. Place a suitable vessel below the discharge tube.
9. The discharge tube must point away from the user at all times.

#### **Warning!**

Only operate the instrument with correctly installed discharge tube. Never push down the piston while the closure cap is mounted!

#### **Note:**

Before using the instrument for the first time, ensure it is rinsed carefully or discard first few samples dispensed.

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1. Hold a suitable vessel below the discharge tube outlet.

**Warning!**  
Reagent may splash!

2. Gently pull up the piston and push down quickly. Repeat this motion in the lower filling area of the cylinder until air bubbles have disappeared from the cylinder.

## 6.3. Setting the volume

Loosen the volume selector screw, set the pointer to the desired volume, and retighten the screw lightly by hand.

## 6.4. Dispensing

1. Verify that the discharge tube is correctly installed and the discharge tube closure cap is removed.
2. Place a suitable vessel below the discharge tube outlet.
3. Gently lift piston in a slow and steady motion, then depress piston slowly and steadily. No bubbles should be visible during dispensing.

**Warning!**  
If the piston moves stiffly or is difficult to move, stop dispensing immediately and clean the instrument (see chap. 7).

4. After dispensing, leave the piston in the down position. Close the discharge tube with the closure cap.

**Attention!**  
Reagent may escape from the discharge tube and the discharge tube closure cap.

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## 7. Cleaning

Instrument will only function safely if cleaned regularly. Be sure to clean the instrument:

1. Immediately, if the piston moves stiffly;
2. Before changing reagents;
3. Before longer downtimes;
4. Before any maintenance and repair work;
5. Before autoclaving.

### Warning!

Dispenser, intake tube, discharge tube closure cap and discharge tube may be filled with reagent. To prevent injury from chemicals, wear eye protection, protective clothing and protective gloves. Always point the discharge tube outlet away from the user and other people.

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### 7.1. Emptying

1. Mount the discharge tube closure cap.
2. Place the instrument with bottle into a suitable basin.
3. Select the maximum volume.
4. Unscrew the dispenser and lift it high enough so that the intake tube (10) is no longer immersed in liquid.
5. Cautiously tap the intake tube against the inside of the bottle so that the reagent runs out.
6. Remove the dispenser from the bottle and mount it on another empty bottle.
7. Remove discharge tube closure cap. Hold the discharge tube outlet over the opening of the first bottle. Empty the dispenser by repeated dispensing movements. For dispensers with a recirculation valve, set the valve switch to "recirculation" and again move the piston repeatedly to empty the recirculation channel.

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## 7.2. Standard cleaning

1. Mount the dispenser on a bottle filled with a suitable cleaning solution.
2. Rinse the dispenser by repeatedly dispensing at maximum volume. For dispensers with a recirculation valve, set the valve switch to "recirculation" and again move the piston repeatedly to rinse the recirculation channel.
3. Empty the dispenser completely, as described in chapter 7.1.
4. Mount the dispenser on a bottle filled with distilled water. Rinse thoroughly, then empty, as described in chapter 7.1.

## 7.3. Intensive cleaning

Intensive cleaning must first be preceded by a standard cleaning. This procedure must be followed if the piston is difficult to move, or if the dispenser is to be autoclaved. For this purpose, the dispenser must be partially disassembled.

### Warning!

Before disassembly, always perform the standard cleaning procedure (see chapter 7.2). To prevent injury from chemicals, always wear eye protection, protective clothing and protective gloves. Avoid splashing reagent.

1. Pull off the intake tube and clean it with a soft bottle brush. Replace as necessary, especially if cracked or worn.
2. Unscrew the piston seat (1) and carefully pull out the piston (4).
3. Clean piston and cylinder with a soft bottle brush and rinse with water.
4. Re-insert the piston carefully so as not to damage the piston seal.
5. Hand tighten the piston seat.
6. Rinse the dispenser thoroughly with distilled water, then empty as described in chapter 7.1.

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## 8. Cleaning/replacing valves

### Warning!

Before disassemble, always perform the standard cleaning procedure (see chapter 7.2). To prevent injury from chemicals, always wear eye protection, protective clothing and protective gloves. Avoid splashing of reagent.

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### 8.1. Cleaning/replacing the intake valve

#### Note:

Valve ball (8) may fall out when valve is removed.

1. Unscrew the intake valve (9) using the valve tool (15).
2. Clean or replace the valve.
3. To reassemble, tighten intake valve (9) by hand and complete using valve tool (15).

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### 8.2. Cleaning/replacing the discharge valve

1. Unscrew the discharge tube (13).
2. Loosen discharge valve securing nut (12).
3. Carefully pull out the discharge valve (11).
4. Clean or replace the valve.
5. To reassemble, push in the valve (11) by hand and tighten discharge valve securing nut (12).
6. Push in the discharge tube and firmly finger-tighten the locking nut.

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### Warning!

Never use force during disassembly and assembly. Verify that all components fit tightly and securely. Subsequent to reassembly, gravimetrically check volume (see chapter 10).

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## 9. Autoclaving

The dispenser can be autoclaved without disassembly (121 °C, 2 bar,  $t_{\max}$  20 min). It is up to the user to verify the effectiveness of the sterilization.

1. Perform intensive cleaning (see chapter 7.3).
2. Pull off the discharge tube closure cap (14) and intake tube (10), and place into the autoclave on a towel, together with the instrument. Avoid contact with hot metallic surfaces.
3. Place instrument in autoclave.

### Attention!

After autoclaving, allow the instrument to slowly cool down to room temperature (approx. 2 hours). Check all parts for deformations, cracks and leaks.

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## 10. Volume Check

The procedure is described in detail in DIN EN ISO 8655/6.

1. Set to nominal volume and dispense distilled water into a weighing vessel.
2. Weigh the dispensed amount of water with an analytical balance.
3. Convert the indicated weight into volume units, taking into account the temperature, density and atmospheric pressure. (The conversion factor Z is published in the standard.)
4. Carry out 10 dispensing procedures.
5. Calculate accuracy and coefficient of variation.

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## 11. Troubleshooting

What to do if ...

1. **Reagent is not aspirated:**
  - Clean intake valve (9); replace if  
necessary (see chapter 8.1).
2. **Cylinder does not fill:**
  - Clean discharge valves (see chapter 8.2).
3. **Visible air bubbles in aspirated liquid:**
  - Check intake valve, make sure properly seated.
  - Check intake tube; replace if necessary.
  - Aspirate reagent more slowly.
4. **Piston movement is stiff:**
  - Clean dispenser thoroughly  
(see chapter 7.3)

## 12. Repair Service

### Attention!

Only cleaned and decontaminated instruments  
can be evaluated/repared.

Clean the instrument thoroughly. If the instrument  
has been contaminated with infectious,  
genetically modified or biologically active or  
radioactive substances, the instrument must be  
decontaminated in addition. Return instrument  
along with a description of the malfunction and  
dispensed liquids.

Return shipment has to be effected at the risk and  
cost of the sender.

## 13. Warranty

We shall not be liable for the consequences of  
improper handling, use, servicing or unauthorized  
repairs of the instrument or the consequences  
of normal wear and tear especially of wearing  
parts such as pistons, seals, valves and the  
breakage of glass as well as the failure to follow  
the instructions of the operating manual. We are  
not liable for damage resulting from any actions  
not described in the operating manual or if non-  
original parts have been used.

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## 14. Technical Data / Ordering Information

### 14.1 Items supplied

Complete with discharge tube, intake tube, 3 adapters and valve tool

- **Nominal volume 2,5 to 10 ml Mounting thread GL 32;**  
Adapters (PP) supplied for bottle threads GL 28, S 40 and GL 45
- **Nominal volume 25 to 100 ml Mounting thread GL 45;**  
Adapters (PP) supplied for bottle threads GL 32, GL 38 and S 40

The dispenser is conformity-certified.

Einstellbereich Volume range	Dosierschritte Dispensing steps	Systematische Messabweichung*, R% Systematic error*, A%	Zufällige Messabweichung*, V% Random error*, CV%
0,5 - 2,5 ml	0,05 ml	≤ ± 0,6	≤ 0,1
1,0 - 5,0 ml	0,10 ml	≤ ± 0,5	≤ 0,1
2,0 - 10,0 ml	0,20 ml	≤ ± 0,5	≤ 0,1
5,0 - 25,0 ml	0,50 ml	≤ ± 0,5	≤ 0,1
10,0 - 50,0 ml	1,00 ml	≤ ± 0,5	≤ 0,1
20,0 - 100,0 ml	2,00 ml	≤ ± 0,5	≤ 0,1

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<b>Calibration conditions</b>	Distilled water, delivered at 20 ± 0.5 °C
<b>Number of testing procedures</b>	10 according to DIN EN ISO 8655/6
(Technical specifications subject to change.)	

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