**Gebrauchsanweisung
Instruction manual
Manual de instrucciones**

**Destillierapparat für Wasser
Distilling apparatus for water
Destilador de agua**

DEST35

- A = Wasseranschluss
Water connection
Conexión de agua
- B = Wasserauslauf
Water outlet
Salida de agua
- C = Destillatauslauf
Distillate drain
Salida de destilado
- D = Entlüftung
Aeration cap
Ventilación
- E = Ablasshahn
Discharge tube
Grifo de descarga

Inhalt

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1. Scope of delivery

- 1 Distilling apparatus for water Cat. no. 7224002, 230V
- 1 Instruction manual

2. General information

The roller mixers are designed in accordance with Safety Class 1. They are produced and tested in accordance with DIN EN 60335-1; VDE 0700-1:2012-10:2012-10. According to these regulations, the units are designed to meet the requirements for safe and correct operations. To ensure the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

2.1. Unpacking the instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport, such as indented sides of the carton, chemical impairment of the carton, water stains or other physical impacts which might have damaged the content, be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage. You may unpack the carton carefully if there is no visible damage of the packing material.

Note:

If the instrument is not damaged and all parts are complete, you may start to operate the device after reading the instruction manual.

Important note:

Read this instruction manual carefully before operating the instrument. Should there be any further questions after reading the manual concerning the setup, operation or warranty, please contact:

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
3. Safety instructions


Important note:


Please store the manual in a place easily accessible for every user.

A EURO AC plug (DIN 49441 CEE 7/VII10/ 16 A 250 V) is normally shipped with our instruments.

When operating the instruments in countries with different AC plug systems, use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation

The instrument is earthed. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!

Caution:

Please comply with all safety and accident-prevention regulations applicable to laboratory work.

Caution:

Instructed users only may operate the instrument.

Caution:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

Caution:

The On/Off switch does not disconnect the instrument from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.

Make sure the unit is standing on a solid and horizontal surface.

4. Intended use

The distilling apparatus is used for mono-distillation of tap-water.

5. Setup and operation

5.1 Setup of the instrument

The attached glass part (packed separately) has to be fitted in the threaded sockets and screwed carefully. Then tubes and mains supply have to be connected.

Caution:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

5.2 Supply

Supply is made by connection A.

First the water runs through the cooling coil, fills the boiling unit and leaks from the water outlet B as soon as reaching the required filling level (automatic level regulation).

The diameter of the drain tube must not be too small since in case of reflux the level regulation does not work faultlessly. An inside diameter of the tube of 12 mm is recommended.

5.3 Heating

The water is brought to boil by means of the heating spiral made of refined steel. A probe at the spiral switches off the heating if there is a defect in the supply of cooling water.

If the temperature falls below the switch-on limit of 110 °C, the heating begins to turn on. This is done automatically as soon as the required water level in the boiling unit is reached again. Already after approx. 3 minutes the water in the boiling unit begins boiling.

5.4 Distillate

The drain of distillate is made by means of the distillate drain C. Please take care that, before putting the appliance into operation, the tubing nipple is pulled out from the distillate as far as possible.

For the first operation discard the first 5 litres of distillate.

To secure the purity of the distillate lead it into a suitable container by means of a PTFE tube. This container is not part of the delivery.

5.5 Draining

To discharge the boiling unit, close the aeration cap D and then close the water outlet B with a suitable stopper or pinchcock. As soon as the heating has switched off, turn on the tap water a bit more for some seconds and then turn it off completely. The content of the boiling unit drains off.

To remove the liquid completely, tilt the appliance to the side and drain off the rest of liquid through the re-opened water outlet B.

Attention: Avoid lateral knocking – risk of glass breakage!

5.6 Cleaning

The lime residues resulting from evaporation of water can be removed with a special lime solvent or 10 % formic acid.

For this purpose the water in the boiling unit has to be drained off (see point 5.5). Insert a glass funnel in the aeration nozzle D, close the water outlet B with a suitable stopper or pinchcock and fill in the lime solvent. Then start re-supply.

To accelerate the procedure of dissolving switch on the heating for a moment or let the solvent take effect overnight.

Attention: Be careful when using a caustic lime solvent.

Draining is to be made as described in point 5.5. Then rinse several times with water (approx. 5 times).

5.7 Maintenance

The appliance does not need any maintenance except the regular removing of the lime residues. Please consider that there is a natural ageing of several elements. In case of any questions, please let us know the type of your equipment and the fabrication number.

6. Dismantling and transport

6.1 Dismantling

1. Switch off the instrument by means of the rocker switch.
2. Disconnect the instrument from the mains.
3. Now you can remove the instrument from the working area.

6.2 Transport

Place the instrument and its parts in its original packaging or in a suitable container to protect it during transport. Close the packaging with adhesive tape. In case you do not use the original packaging, please mark the box with the following notes:

- Glass symbol (Fragile! Handle with care!)
- Umbrella (Keep dry!)
- Content (List of content)

7. Disposal

Damaged and/or waste electric or electronic equipment have to be disposed of at the designated Recycling Depot. Also the packing material should be disposed of environment-friendly (material separation)

8. Warranty and Liability

Paul Marienfeld GmbH & Co. KG provides a guarantee for a period of 24 months from date of delivery for the instrument purchased by you. This warranty applies only to defects in material or workmanship. There is no warranty for defects or malfunctions caused by neglect, improper use or improper service.

In case of defects in material or workmanship, the instrument will be repaired or replaced. In case of replacement, the warranty for the new unit will persist for the remaining period of the primary warranty period. Charges for returning the instrument have to be borne by the customer.

Any claims regarding guarantee and liability are limited to professional and proper use of the instrument and observance of the general safety instructions and this instruction manual.

Paul Marienfeld GmbH & Co. KG is not liable for any material or personal damages caused by improper use and does not accept liability for consequential damages. The warranty expires in case of removal or alteration of the identification label and/or the serial number of the instrument.

9. Technical Data

Mains:	230 V/AC
Power consumption:	2500 W
Output of distillate:	approx. 3.5 l/h
Required consumption of cooling water:	approx. 0.75 l/min
Outlet temperature of cooling water while inlet temperature is 18 °C:	approx. 60 °C
Conductivity of distillate at 20 °C:	0.85 µS
Material of tube heating:	refined steel 1.4541
Material of glass part:	borosilicate glass 3.3
Total weight:	approx. 4 kg

10. Declaration of conformity

We hereby declare that this product

Distilling apparatus for water DEST35

is in conformity with the following directives:

EMV directive 89/336/EEC amended by 91/263/EEC

EMV directive 72/23/EEC

and is in accordance with the following standards or standardized documents:

EN 292-1/11.91

EN 292-2/11.91

Electromagnetic compatibility:

EN 50081-1/3.93

EN 50082-1/11.97

In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.

Paul Marienfeld GmbH & Co. KG
Harry Marienfeld, President¹

¹ The above declaration has been generated automatically and is valid without signature

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