**Gebrauchsanweisung
Instruction manual
Manual de instrucciones**

**Laborrührwerk R14
Laboratory stirrer R14
Agitador de laboratorio R14**

Inhalt

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Content

1. Scope of delivery	13
2. General information	13-14
2.1 Unpacking the instrument	13, 14
3. Safety instructions	14-16
3.1 Description of safety symbols in the instruction manual	14
3.2 Description of safety symbols on the instrument	15
3.3 General safety instructions	15, 16
4. Intended use	16
5. Setup and operation	17
6. Cleaning and maintenance	17, 18
7. Dismantling	18
8. Transport and storage	18
9. Disposal	18
10. Warranty and liability	19
11. Technical data	19

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optionally available:

- Cat. No. 6641301 Blade stirrer
- Cat. No. 6647402 Support plate for use with support rod
- Cat. No. 6647601 Support rod for use with support plate
- Cat. No. 6628000 Bosshead to fasten the stirrer to a stand

2. General information

Our laboratory stirrers are designed according to the protection class 1. They were manufactured and tested according to DIN EN 61010. According to these regulations, the units are designed to meet the requirements for safe and correct operations. To ensure the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

2.1. Unpacking the instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

Note:

If the instrument is not damaged and all parts are complete, you may start to operate the device after reading the instruction manual.

Important note:

Read this instruction manual carefully before operating the instrument. Please store the manual in a place easily accessible for every user. Should there be any further questions after reading the manual concerning the setup, operation or warranty, please contact:

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A EURO AC plug (DIN 49441 CEE 7/VII/10/ 16 A 250 V) is normally shipped with our instruments.

When operating the instruments in countries with different AC plug systems, use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation. Ensure that your local supply voltage matches the indication on the instrument's specification plate!

The instrument is earthed. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!

3. Safety instructions

3.1. Description of safety symbols in the instruction manual

Warning against general danger: This symbol indicates that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user.

This symbol distinguishes notices which have to be strictly observed by the user. Disregard of warnings may result in the impairment of serviceability as well as in physical harm to the user.

This symbol distinguishes notices which have to be strictly observed by the user to ensure safe operation of the unit.

Caution
Danger for fire or explosion!

Note regarding repair or maintenance

Note regarding mains cable connection

3.2. Description of safety symbols on the instrument

Warning against general danger: This symbol indicates that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user.

3.3. General safety instructions

Important note:

It is under the responsibility of the user to find out whether the device is suitable for his application. If in doubt, clarify this with your dealer or with Paul Marienfeld GmbH & Co. KG directly.

Caution:

Please comply with all safety and accident-prevention regulations applicable to laboratory work.

Caution:

Please do not use this instrument in insecure environment especially not in explosive ambience. Danger of life!

Caution:

Instructed users only may operate the instrument.

Caution:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

Note:

Make sure that all safety information on the instrument is clearly visible during operation.

Attention:

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

The mains plug may only be inserted in a socket with earthing contact. The protective effect must not be nullified by an extension cable without protective conductor. Any interruption of the protective conductor inside or outside the device or disconnection of the protective conductor connection may cause the device to become dangerous. Intentional interruption is not permitted.

To ensure safe operation fasten the instrument on a support stand.

The On/Off switch does not disconnect the instrument from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.

Protect the vessel containing the liquid to be stirred against rotation, vibration and movement by means of suitable clamping devices.

Before changing the stirrer blade first switch off the instrument and disconnect the instrument from the mains supply entirely.

Do not start the motor until the properly mounted stirring blade is immersed deep enough in the liquid.

Make sure that the liquid cannot be splashed out of the vessel by using a suitable container which is high enough to make sure the walls protrude adequately over the liquid level. A suitable lid is useful as a splash guard.

Wear protective clothing and gloves in case you expect splashing of the liquid.

Make sure the chuck is not immersed into the liquid to be mixed.

4. Intended use

These laboratory stirrers are mainly to be used to stir liquids of low to medium viscosity. They are designed for use in chemical and biological laboratories or industrial enterprises, universities and pharmacies.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.

5. Setup and operation

Caution:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

Caution

Please do not use this instrument in insecure environment especially not in explosive ambience. Danger of life!

Caution:

Disconnect the instrument from the mains supply as a protection against unauthorised or accidental operation whenever you stop the work.


Caution:

Always start stirring with the lowest speed and then rev up according to the required intensity.

The laboratory stirrer R14 is a convenient, robust overhead stirrer for use with 4 mm diameter stirring rods.

1. Mounting the support rod:
Screw the support rod clockwise into the thread bore at the back of the stirrer. Use a suitable 6 mm tool to firmly tighten the support rod.
2. Mounting the stirring tool:
 - a. Unfasten the keyless chuck (B)
 - b. Move the piston of the stirring tool (C) into the chuck
 - c. Torque down the chuck (B) when the piston is in the intended position.
3. Insert the plug into the mains outlet.
4. Place the stirring rod into the media to be stirred.
5. Switch on the instrument with the speed control knob and choose the intended speed.
The stirring motor is designed for continuous operation. It has an electronically controlled speed of 100 to 2000 rpm.

6. Cleaning and maintenance

The outer casing of the instrument consists of lacquered or anodised aluminium and a chemically resistant splash-proof membrane key pad. To clean the instrument use a damp cloth and a mild soap solution. Do not use steel wool or any similar aggressive cleaning agents to clean the instrument. The surface of the instrument could be damaged.

**Attention:**

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

The stirrer is maintenance free and only subject to a natural ageing process and the statistical failure rate of the components. All bearings of the stirring motor are sufficiently lubricated so that with careful handling no maintenance is necessary. Should there be any problems, please send the instrument to Marienfeld or your agent.

**Advice for service technicians only:**

The device is 2-phase protected. The electronic fuses are located on the PCB. Disconnect the instrument from the mains supply before opening the instrument.

7. Dismantling

Switch of the stirrer and disconnect from the mains. Remove any glass beakers and any other equipment around the overhead stirrer.

Now unfasten the keyless chuck and remove the stirrer blade. After that remove the stirrer from the support and remove the support rod.

8. Transport and storage

Please use, if possible, the original packaging. Protect the instrument by means of bubble wrap against mechanical shocks from the outside. In case you do not use the original packaging, please mark the box with the following notes:

- Glass symbol (Fragile! Handle with care!)
- Umbrella (Keep dry!)
- Content (List of content)

For proper storage please observe the following ambient conditions:

Max. temperature: +5 °C bis +40 °C
 Max. humidity: 80 %

9. Disposal

Damaged and/or disposed electric or electronic equipment has to be disposed of at the designated recycling depot. Also the packing material should be disposed of environment-friendly (material separation).

10. Guarantee and liability

Paul Marienfeld GmbH & Co. KG provides a guarantee for a period of 24 months from date of delivery for the instrument purchased by you. This warranty applies only to defects in material or workmanship. There is no warranty for defects or malfunctions caused by neglect, improper use or improper service.

In case of defects in material or workmanship, the instrument will be repaired or replaced. In case of replacement, the warranty for the new unit will persist for the remaining period of the primary warranty period. Charges for returning the instrument have to be borne by the customer.

Any claims regarding guarantee and liability are limited to professional and proper use of the instrument and observance of the general safety instructions and this instruction manual.

Paul Marienfeld GmbH & Co. KG is not liable for any material or personal damages caused by improper use and does not accept liability for consequential damages. The warranty expires in case of removal or alteration of the identification label and/or the serial number of the instrument.

11. Technical data

Mains:	115 V resp. 230 V, 50-60 Hz (see rear panel)
Motor:	25 Watt, electronically controlled speed
Speed:	100 to 2000 rpm
Weight:	approx. 1.2 kg
Dimensions:	66 x 120 x 64 mm (L x W x H)
Mounting support:	190,5 x 16 mm

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