## Inhalt

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## Content

1. Scope of delivery	17
2. General information	17
2.1 Unpacking the instrument	17
2.2 Safety instructions	17, 18
3. Intended use and description	18
3.1 Intended use	18
3.2 Description	18, 19
4. Setup	19
5. Operation	20-22
6. Installation of the software	22, 23
7. Application of the software	23
7.1 Starting the application	23
7.2 Reading counter data and saving patient information	23
7.3 Loading patient information	24
7.4 Backup of patient information database	24
8. Additional functions	24
8.1 Only one Dataset	24
8.2 Multi Count	24
8.3 Send Data to Clipboard	25
9. Troubleshooting	25
10. Cleaning and maintenance	25
11. Dismantling and disposal	26
12. Transport and storage	26
12.1 Packing and transport	26
12.2 Storage	26
13. Warranty and liability	26, 27
14. Technical data	27

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## 1. Scope of delivery

- 1 Cellcounter 2001  
Cat. no. 7140003 (230V/50-60Hz) or Cat. no. 7110003 (115V/50-60Hz)
- 1 Instruction manual
- 1 Windows Software Counter Diff 15 PC Manager

## 2. General information

The instrument is designed in accordance with safety class II. It is built and tested in accordance with DIN EN 60335-1; VDE 0700-1:2012-10:2012-10, safety regulations for electrical instruments. It left our warehouse working perfectly well and safely. To maintain the proper safety and operational function of the instrument, the user must respect the instructions and safety guidelines of this manual.

The instrument comes with a SCHUKO plug (DIN 49441, 16 A, 250 V). When operating the instrument in countries with different AC plug systems, use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation.

### 2.1 Unpacking the instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport, such as indented sides of the carton, chemical impairment of the carton, water stains or other physical impacts which might have damaged the content, be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage. You may unpack the carton carefully if there is no visible damage of the packing material.

### 2.2 Safety instructions

#### Note:

If the instrument is not damaged and all parts are complete, you may start to operate the device after reading the instruction manual.

#### Important note:

Read this instruction manual carefully and make sure everything is understood before operating the instrument. Should there be any further questions after reading the manual concerning the setup, operation or warranty, please contact:

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**Attention:**

Please comply with all safety and accident-prevention regulations applicable to laboratory work.

**Caution:**

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

**Caution:**

The On/Off switch does not disconnect the instrument from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.

**Attention:**

Do not open the instrument! Repairs are only to be carried out by trained service technicians. Any work on the electronics of the unit should only be done by knowledgeable and trained personnel.  
In the case of unauthorized operation as well as negligent or wilful damage, the warranty expires.


### 3. Intended use and description

#### 3.1 Intended use

The Cellcounter 2001 is an electronic memory counter for statistical evaluation. It is especially designed for the counting of blood cells but is also appropriate for counting any other particles.

#### 3.2 Description

The Cellcounter 2001 is an electronic count instrument with 12 stores, counting keys (4) with light indicators, four defined programme keys (3) and a free selection programme, LED display (2) and a rocker switch (1).

The upper two key rows consist of programming and operating keys only. The lower three key rows consist entirely of counting keys.

The instrument has 4 programme keys:

**LEUKO**  
**MYELO / RBC**  
**RETI**  
**FREE** (free selection programme)

and 6 operational keys:

**SET**  
**GO ON + COUNT**  
**%**  
**REAL**  
**CLEAR**  
**ERROR CORR**

#### 4. Setup

Please observe the ambient conditions indicated in the technical data (page 20). Make sure that your local supply voltage corresponds to the specifications indicated on the instrument.

## 5. Operation

The counter has a self-testing facility that starts automatically when it is switched on. After approx. 2 seconds sounds an acoustic signal, a "0" is displayed and the programme keys LEUKO and CLEAR light up. The Cellcounter 2001 starts automatically with the programme **LEUKO** because this is the most common programme.

For starting to count with **LEUKO** press **SET** repeatedly in order to choose the desired sum of the counting process. It can be set individually from 100, 200, 300 ... up to 1000 → "000"). Only for **Div/Reti** and **Div/Ery** the sum is 1000 by default.

In order to correct possible errors press **CLEAR** for at least 3 seconds. All information previously entered will be cancelled, the counter returns to its start-position and can be programmed once more. Pressing the corresponding key you can choose any of the 4 counting programmes. The total sum can also be programmed pressing **SET**. Having chosen the programme and the counting sum, counting can be started by pressing **GO ON + COUNT**. Now the program can only be changed using **CLEAR** as mentioned above.

For every counted cell the corresponding key has to be pressed once and the corresponding counter counts one impulse. While keeping pressed the key the display shows the content of this counter that means the number of the respective cells. When releasing the key the total sum of all counting impulses or the number of all counted cells. This way the counting can be checked at any time before reaching the predetermined sum of the counting process

During the counting process the percentage can also be displayed at any time by pressing **%**. All keys that have been used for counting will be flashing. The key of which the counting content is displayed at the moment is illuminated. The display shows the percentage of this cell type in relation to the total count. Each of the blinking keys can be chosen in order to display its percentage. Thereafter the light of the checked key will switch off. If you press a key that was not used for counting before, the display shows "--". Counting can be continued with **GO ON + COUNT**.

In case of an error the last entered impulse can be deleted using **ERROR CORR**. Pressing **ERROR CORR** the value of the last pressed counting key is displayed reduced by the last entered value. During this process the respective counting key will be illuminated. Releasing **ERROR CORR** the new counting total will be displayed, reduced by the last entered value of the mentioned register key.

Only the last entered counting impulse can be cancelled.

Reaching the pre-defined sum, an acoustic signal sounds and all entry keys are blocked. The keys **REAL** and **%** can be used as previously described. Further counting may only be carried out after pressing the **GO ON + COUNT** key, otherwise a warning signal sounds.

The 4 programme keys:

**LEUKO** Up to 12 different cell counting keys can be included in the percentage proportion for leukocyte differentiation. The total sum can be selected between 100 and 1000 in stages of 100. Furthermore 3 single count keys (Normo, Tox Gr. and Div/Ery) are available but those cannot be included in the percentage calculation.

**MYELO/RBC** Up to 15 different cell counting keys can be included in the percentage proportion. The total sum can be selected between 100 and 1000 in stages of 100.

**RETI** Only the Div/Reti and Div/Ery keys are available in this programme. The counting steps for the erythrocytes can be set to 1/5/10/15/20/30/40 and 50 before starting the counting process. The total sum is 1000 by default. The reticulocytes are counted in steps of 1. Div/Ery is always 100%. Div/Reti cannot be displayed in percentage – the display shows “- - -”.


**FREE** Within this optional programme all cell counting keys may be either included in or excluded from the percentage calculation. Keys that should not be included in the percentage calculation have to be pressed once before starting to count. Those keys are illuminated afterwards. Pressing them again switches off the light and includes the counting key into the percentage calculation. After determination of the keys the total sum has to be chosen with **SET**, as previously described. Pressing **SET** repeatedly increases the total sum each time for 100 counting impulses till max. 1000. The counting process can be started with **GO ON + COUNT**. Pressing a key which is not included into the percentage calculation a slight acoustic signal sounds. Checking the percentage for those keys with %, “- - -” is displayed while for the other keys the percentage is displayed as usual.

The 6 function keys are:

- SET** This key presets the sum of the count.
- GO ON + COUNT** This key has to be pressed in order to start counting, go on counting after checking the percentage values and in order to go on counting after reaching the predetermined total counting cum.
- %** Displays percentage of the content of the cell counting key in relation to the whole counting.
- REAL** Displays the absolute content of a cell counting key.
- CLEAR** Resets the counter to "0". In order to avoid equivocated pressing of this key, it has to be pressed at least 3 seconds. Now the counter can be programmed once more.
- ERROR CORR** Deletes the last entered counting impulse.

### 6. Installation of the software

All rights at the software programme Counter Diff 15 PC Manager belong to Ingenieurbüro Zipperer. It may be used only with the Cellcounter

The software Counter Diff 15 PC Manager requires Microsoft Windows XP SL3 or a higher version. Before installing the software, make sure that you have also installed the latest service packs for Windows.

1. First, uninstall any earlier versions of the Counter Diff 15 PC Manager if applicable.
2. Insert the CD Counter Diff 15 PC Manager into your CD drive.
3. Open Windows Explorer and select the CD „Counter“.
4. To run the installation programme go to the directory „Counter Diff 15 PC Manager“ on the CD, open the file „setup.exe“ and follow the on-screen instructions.
5. This programme requires Microsoft .NET Framework 2.0. If not already installed, it will be done automatically.

## 7. Application of the software

### 7.1 Starting the application

1. First, connect the Cellcounter to the serial interface of the PC, then to the mains.
2. Switch on the Cellcounter.
3. Double-click the icon Counter Diff 15 PC Manager on the desktop of the computer.

### 7.2 Read Out Counter Data and save patient information

1. If you have connected the counter to any other COM port than COM port 1 then go to the menu Settings and then to menu COM port to set the appropriate COM port.
2. To change a button name, please click on the button name in the software and alter as required.
3. Start counting with the Cellcounter.
4. Having finished counting, click the **Read Out Counter Data** button in the software to upload the data to the manager. If this does not work, please refer to the trouble shooting section.
5. Next clicking the **Save Data** button will open a window to enter patient information (last- and first name are required). When finished, click the **Next** button to confirm saving all entered data. The counter data and the patient information will be saved to disk.
6. To clear the counter for reading the next data set, press the **Clear** button of the Cellcounter and hold for three seconds. Now continue from step 3.

### 7.3 Loading patient information

If you want to see saved data of a filed patient again, please click the **Search Data** button. A new window appears where you can search for the patient. To see the corresponding counter data, click the **OK** button. The search window will close and the relevant counter data is displayed.

The data can be printed at this point if required. Select the **File** menu and then **Print**.

### 7.4 Making a backup of patient information database

To backup all saved patient data, go to the menu **File** and click **Open Database Directory**. After that you will see a database file which can be saved with a programme of your choice.

## 8. Additional functions

### 8.1 Only one Dataset

**Caution:**

If patient information is already saved in the programme, make a backup of the database before you continue (see 7.4)

If you want to use the saved data of a filed patient in another software (laboratory database, Excel), many of them need only one dataset. This dataset could be loaded for further processing. Therefore, activate the function **Only One Dataset** in menu **Settings**.

**Caution:**

After activating the function **Only one Dataset**, all saved patient information will be permanently deleted.

To see the location of the database you have to select the **File** menu and then **Open Database Directory**.

### 8.2 Multi Count

If you want to compare different counts quickly, activate the function **Multi Count** in the **Settings** menu. The buttons will change and you can navigate between the different counts and to complete them. The function **End with....** in the menu **Settings** offers the choice between printing or saving the data after having activated the button **End**.

### 8.3 Send Data to Clipboard

To easily transfer counter data into another program (i.e. spreadsheet), the Windows clipboard can be used. Activate the function **Send Data to Clipboard** in the menu **Settings**. After reading out counter data, these data will be pasted automatically into the Windows clipboard. Afterwards, counter data can be inserted into the programme of your choice.

### 9. Troubleshooting

If you are uncertain which COM port you use, please click the **Read Out Counter Data** button to check whether the communication works. If not, repeat the setting with COM port 2, COM port 3, etc. from the **Settings** menu.

### 10. Cleaning and maintenance

The case is made of acid-resistant material and easy to clean. Surface and operating elements may be cleaned with a mild dishwashing detergent (water and a standard dishwashing detergent) and a soft, lint-free **moist** cloth. Under no circumstances the surface should be cleaned with acids, chlorinated solvents or saline solutions. Do not rub or scratch to avoid damages of the inscriptions and the surface. If you use any other cleaning method, please make sure that the intended method does not cause any damage to the instrument.

**Attention:**

Carefully prevent that no liquid penetrates into the spaces between the keys. Do not rinse the instrument.

## 11. Dismantling and disposal

1. Switch off the instrument by means of the rocker switch.
2. Disconnect the instrument from the mains.
3. Now you can remove the instrument from the working area.

Damaged and/or disposed electric or electronic equipment has to be disposed of at the designated recycling depot.

Also the packing material should be disposed of environment-friendly (material separation).

## 12. Transport and storage

### 12.1 Packing/Transport

Please use, if possible, the original packaging. Protect the instrument by means of bubble wrap against mechanical shocks from the outside. In case you do not use the original packaging, please mark the box with the following notes:

- Glass symbol (Fragile! Handle with care!)
- Umbrella (Keep dry!)
- Content (List of content)

### 12.2 Storage

For proper storage please observe the following conditions of surrounding:

max. ambient temperature: 5-40 °C

max. ambient humidity: 80 %

## 13. Warranty and Liability:

### Warranty:

Paul Marienfeld GmbH & Co. KG provides a guarantee for a period of 24 months from date of delivery for the instrument purchased by you. This warranty applies only to defects in material or workmanship. There is no warranty for defects or malfunctions caused by neglect, improper use or improper service.

In case of defects in material or workmanship, the instrument will be repaired or replaced. In case of replacement, the warranty for the new unit will persist for the remaining period of primary warranty period. Charges for returning the instrument must be borne by the customer.

Any claims regarding guarantee and liability are limited to professional and proper use of the instrument and observance of the general safety instructions and this instruction manual.

### Liability:

Paul Marienfeld GmbH & Co. KG is not liable for any material or personal damages caused by improper use and does not accept liability for consequential damages. The warranty expires in case of removal or alteration of the identification label and/or the serial number of the instrument.

#### Attention:

In case of malfunction do not try to carry out any repair works. The instrument does not contain any part which may be serviced or maintained by the user. In case of any attempt by the user to repair the unit, the warranty expires!

Do not open the instrument. Any work on the electronics of the unit should be done by knowledgeable and trained personnel only!

## 14. Technical data

Power supply:	115 Volt or 230 Volt (see back panel)
Power input:	6 W
Nominal value of fuse:	250 mA
Permissible ambient temperature:	5-40 °C
Permissible relative humidity:	80 %
Admissible operation:	100 %
Protection class DIN EN 61140:	II
Protection class DIN 40050:	IP 40
Dimensions (LxWxH):	210x190x25 mm
Weight:	0.89 kg

## 1. Alcance de suministro

1 Cellcounter 2001

Cat. n° 7140003 (230V/50-60Hz) o Cat. n° 7110003 (115V/50-60Hz)

1 Manual de instrucciones

1 Windows Software Counter Diff 15 PC Manager

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