**Gebrauchsanweisung
Instruction manual
Manual de instrucciones**

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Content

1. Scope of delivery	15
2. General information	15-16
2.1 Unpacking the instrument	15, 16
3. Safety instructions	16-18
3.1 Description of safety symbols in the instruction manual	16
3.2 Description of safety symbols on the instrument	17
3.3 General safety instructions	17, 18
4. Intended use	19
5. Mounting and setup	19-20
5.1 Mounting the support rod	19
5.2 Mounting the stirrer on a stand	19
5.3 Mounting the stirring paddle	20
5.4 Mounting the shaft protector	20
6. Operation	21
7. Cleaning and maintenance	21, 22
8. Dismantling	22
9. Transport and storage	22, 23
10. Warranty and liability	23
11. Technical data	24

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optionally available:

- Cat. no. 6641302 Blade stirrer
- Cat. no. 6641303 Propeller stirrer
- Cat. no. 6647402 Support plate for use with support rod
- Cat. no. 6647601 Support rod for use with support plate
- Cat. no. 6628000 Bosshead to fasten the stirrer to a stand
- Cat. no. 7205227 Shaft protector

2. General information

Our laboratory stirrers are designed according to the protection class 1. They were manufactured and tested according to DIN EN 61010. According to these regulations, the units are designed to meet the requirements for safe and correct operations. To ensure the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

2.1. Unpacking the instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

Note:

If the instrument is not damaged and all parts are complete, you may start to operate the device after reading the instruction manual.

Important note:

Read this instruction manual carefully before operating the instrument. Please store the manual in a place easily accessible for every user. Should there be any further questions after reading the manual concerning the setup, operation or warranty, please contact:

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A SCHUKO plug (DIN 49441, 16 A, 250 V) is normally shipped with our instruments.

When operating the instruments in countries with different AC plug systems, use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation. Ensure that your local supply voltage matches the indication on the instrument's specification plate!

The instrument is earthed. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!

3. Safety instructions

3.1. Description of safety symbols in the instruction manual

Warning against general danger: This symbol indicates that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user.

This symbol distinguishes notices which have to be strictly observed by the user. Disregard of warnings may result in the impairment of serviceability as well as in physical harm to the user.

This symbol distinguishes notices which have to be strictly observed by the user to ensure safe operation of the unit.

Caution
Danger for fire or explosion!

Note regarding repair or maintenance

Note regarding mains cable connection

3.2. Description of safety symbols on the instrument

Warning against general danger: This symbol indicates that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user.

3.3. General safety instructions

Important note:

It is under the responsibility of the user to find out whether the device is suitable for his application. If in doubt, clarify this with your dealer or with Paul Marienfeld GmbH & Co. KG directly.

Caution:

Please comply with all safety and accident-prevention regulations applicable to laboratory work.

Caution:

Please do not use this instrument in insecure environment especially not in explosive ambience. Danger of life!

Caution:

Instructed users only may operate the instrument.

Caution:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

Note:

Make sure that all safety information on the instrument is clearly visible during operation.

Attention:

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

The mains plug may only be inserted in a socket with earthing contact. The protective effect must not be nullified by an extension cable without protective conductor. Any interruption of the protective conductor inside or outside the device or disconnection of the protective conductor connection may cause the device to become dangerous. Intentional interruption is not permitted.

To ensure safe operation fasten the instrument on a support stand.

The On/Off switch does not disconnect the instrument from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.

Protect the vessel containing the liquid to be stirred against rotation, vibration and movement by means of suitable clamping devices.

Before changing the stirrer blade first switch off the instrument and disconnect the instrument from the mains supply entirely.

Do not start the motor until the properly mounted stirring blade is immersed deep enough in the liquid.

Make sure that the liquid cannot be splashed out of the vessel by using a suitable container which is high enough to make sure the walls protrude adequately over the liquid level. A suitable lid is useful as a splash guard.

Wear protective clothing and gloves in case you expect splashing of the liquid.

Make sure the chuck is not immersed into the liquid to be mixed.


4. Intended use

The microprocessor controlled overhead stirrers R18 and R50 resp. R50 D are designed for qualified use in laboratories. They are mainly used to stir liquids from aqueous to medium viscosity. For correct use the instrument must be fixed to an appropriate stand.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.


5. Mounting and setup

5.1 Mounting of the support rod

- Fasten the support rod (A) by screwing the two M5 screws (B) into the provided threaded holes at the rear side of the panel.
- Tighten the screws thoroughly using a suitable hexagon key.

5.2 Mounting the stirrer on a stand

- Fasten the bosshead on the rod of the stand (C) at one end.
- Insert the support rod at the other end of the bosshead (D)
- Tighten the thumbscrew (E) until the support rod is securely fastened


5.3 Mounting the stirring paddle

The stirrers R50 and R50 D are equipped with a hollow shaft of 10.5 mm diameter. The chuck is chromium-plated for prevention of corrosion. It takes any stirring rod up to 10 mm diameter. The support rod (A) included in the delivery is used to attach the motor to a stand. The stirrer R18 does not have a hollow shaft but also takes stirrers with a shaft diameter of up to 10 mm.

- Insert the stirring paddle into the chuck and firmly tighten it by hand in the desired position.
- The immersion depth of the stirring tool may be adjusted by loosening the chuck and axially moving the stirring paddle.

5.4 Mounting the shaft protector (optional accessory)

We recommend to use the shaft protector for secure working with the stirrer (optional accessory Marienfeld no. 7205227).

Attach the fastening rods (H) to the stand rod in front of the shaft protector using for example the Jumbo clamp (Marienfeld no. 7205228).

Position the bores (K) of the protector next to the screws (J), insert the screws (J) into the bores and turn the protector to the right as far as it will go.

The screw (L) is used to adjust the overall height of the protector.

Check that the stirring shaft protector is held in position securely prior to each use and also at regular intervals.

The position of the protector must be adjusted only when the stirrer is turned off and the power supply is disconnected.

6. Operation

Caution:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

Caution

Please do not use this instrument in insecure environment especially not in explosive ambience. Danger of life!

Caution:

Disconnect the instrument from the mains supply as a protection against unauthorised or accidental operation whenever you stop the work.

Caution:

Always start stirring with the lowest speed and then rev up according to the required intensity.

There are the following operating elements:

- Toggle switch I/O (on/off)
- R18 and R50: rotary knob with scale for adjustment of the speed
- R50 D: rotary knob for adjustment of the speed, speed will be shown in the digital display

Select the suitable stirring tool for the intended purpose, insert it into the chuck and fasten tight. Set the speed by means of the rotary knob to the lowest possible speed. Switch on the instrument and adjust the speed to the desired value.

Should an overload occur from overheating during continual operations or of a too hot ambient, the red LED will show overload and the electronic will reduce the power intake. The motor will switch off if continuously over heated. After a cool down period and switching the unit off and on again, the motor will operate as normal.

7. Cleaning and maintenance

The outer casing of the instrument consists of lacquered or anodised aluminium and a chemically resistant splash-proof membrane key pad. To clean the instrument, use a damp cloth and a mild soap solution. Do not use steel wool or any similar aggressive cleaning agents to clean the instrument. The surface of the instrument could be damaged.

Attention:

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

Apart from the motor carbon brushes the stirrer is maintenance free and only subject to a natural ageing process and the statistical failure rate of the components. All bearings of the stirring motor are sufficiently lubricated so that with careful handling no maintenance is necessary. Should there be any problems, please send the instrument to Marienfeld or your agent.

Advice for service technicians only:

The device is 2-phase protected. The electronic fuses are located on the PCB. Disconnect the instrument from the mains supply before opening the instrument.

The carbon brushes in the motor are subject to natural abrasion. They can be changed as described following herewith:

1. Switch off the instrument.
2. Unplug the AC power.
3. Remove the black cover at the side of the motor.
4. Bent up the metal clip and take out the spring.
5. Use a small screwdriver to remove the clip from the carbon and take out the carbon.
6. Replace the carbon, attach the spring, bent back the clip and attach the black cover.

In case that the tooth belt needs to be changed, please follow this instruction:

1. Switch off the instrument.
2. Unplug the AC power.
3. Remove the hard plastic undercover.
4. Remove the tooth belt from the large cogwheel.
5. Pull the new tooth belt onto the small cogwheel first and then onto the large cogwheel.
6. Now turn the large cogwheel slightly until the belt is in the proper position.

8. Dismantling

Switch of the stirrer and disconnect from the mains. Remove any glass beakers and any other equipment around the overhead stirrer.

Now unfasten the keyless chuck and remove the stirrer blade. After that remove the stirrer from the support and remove the support rod.

9. Transport and storage

Please use, if possible, the original packaging. Protect the instrument by means of bubble wrap against mechanical shocks from the outside. In case you do not use the original packaging, please mark the box with the following notes:

- Glass symbol (Fragile! Handle with care!)
- Umbrella (Keep dry!)
- Content (List of content)

For proper storage please observe the following ambient conditions:

Max. temperature: +5 °C bis +40 °C
 Max. humidity: 80 %

10. Disposal

Damaged and/or disposed electric or electronic equipment has to be disposed of at the designated recycling depot. Also the packing material should be disposed of environment-friendly (material separation).

11. Guarantee and liability

Paul Marienfeld GmbH & Co. KG provides a guarantee for a period of 24 months from date of delivery for the instrument purchased by you. This warranty applies only to defects in material or workmanship. There is no warranty for defects or malfunctions caused by neglect, improper use or improper service.

In case of defects in material or workmanship, the instrument will be repaired or replaced. In case of replacement, the warranty for the new unit will persist for the remaining period of the primary warranty period. Charges for returning the instrument have to be borne by the customer.

Any claims regarding guarantee and liability are limited to professional and proper use of the instrument and observance of the general safety instructions and this instruction manual.

Paul Marienfeld GmbH & Co. KG is not liable for any material or personal damages caused by improper use and does not accept liability for consequential damages. The warranty expires in case of removal or alteration of the identification label and/or the serial number of the instrument.

12. Technical data

Mains:	115 V resp. 230 V, 50-60 Hz (see rear panel)
Motor:	Permanently generated direct current motor
Speed:	110 to 2000 rpm (R18) 50 to 1500 rmp (R50 and R50 D)
Power consumption:	120 Watt, short term overload 30 %
Torque over entire speed range:	18 Ncm (R18) 55 Ncm (R50 and R50 D)
Capacity of chuck:	0.5 to 10.5 mm
Viscosity of media:	R18 – aqueous fluids R50 and R50 D medium to wallpaper paste
Weight:	approx. 2.1 kg
Dimensions:	65 x 108 (plus 55 mm chuck) x 120 mm (L x W x H)
Diameter of support rod:	16 mm

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