
Gebrauchsanweisung
Instruction manual
Manual de instrucciones

Destillierapparat für Wasser
Distilling apparatus for water
Destilador de agua
DP 4000

Rückansicht
 Rear view
 Vista trasera

Vorderansicht
 Front view
 Vista frontal

Inhalt

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Content

1. Scope of delivery	19
2. General information	19
2.1 Unpacking	19
3. Safety instructions	20, 21
4. Intended use	21
5. Setup and operation	21
5.1 Setup of the instrument	21
5.2 Installing a reservoir bottle	21, 22
5.3 Protective cap	22
5.4 Connecting the tubing	23
5.4.1 Connection of the cooling water	23
5.4.2 Connection of the water drainage	23
5.4.3 Connection of the outlet valve	24
5.5 Preparation of operation	24
5.6 Adjusting the water supply	24
5.7 Interruption of water supply	25
5.8 To switch off	25
5.9 Cleaning	25-29
5.9.1 Cleaning the evaporator	26, 27
5.9.2 Cleaning of the cooling coil	28, 29
6. Dismantling and transport	29
6.1 Disassembly	29
6.2 Transport	29
6.3 Reassembly	30
7. Disposal	30
8. Warranty and liability	30
9. Technical data	31
10. Declaration of conformity	32

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1. Scope of delivery

- 1 Distilling apparatus for water Cat. no. 7224003, 230V
- 1 Instruction manual
- 1 Allen key, hexagonal, 3 mm
- 1 set of tubing:

Material	dia. ext.	dia. int.	wall thickness	overall length
PVC	19 mm	14 mm	2.5 mm	3.5 m
Silikon	16 mm	12 mm	2.0 mm	0.4 m

2. General information

The apparatus is designed in accordance with Safety Class 1. It is produced and tested in accordance with DIN EN 60335-1; VDE 0700-1:2012-10:2012-10. According to these regulations, the unit is designed to meet the requirements for safe and correct operations. To ensure the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

2.1. Unpacking the instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport, such as indented sides of the carton, chemical impairment of the carton, water stains or other physical impacts which might have damaged the content, be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

You may unpack the carton carefully if there is no visible damage of the packing material. Carefully remove the upper part of the interior cover made from polyurethane foam and lift the apparatus out of the carton. Now check the device immediately for glass breakage and other damages that may have occurred during transport.

Note:

If the instrument is not damaged and all parts are complete, you may start to operate the device after reading the instruction manual.

Important note:

Read this instruction manual carefully before operating the instrument. Should there be any further questions after reading the manual concerning the setup, operation or warranty, please contact:

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
3. Safety instructions


Important note:

Please read and follow this instruction manual carefully step by step. It should be available at all times, especially to anyone who wants to use this device. Therefore, this instruction manual should be kept near the device.

The area in which the device is installed must comply with the relevant electrical regulations.
The instrument is earthed. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!

The electrical connection must be protected by a series-connected fuse with at least 15 A (max. 16 A).
The connection to the power supply and the water connection must be made in such a way that there is no danger from dripping water.

Caution:

Please comply with all safety and accident-prevention regulations applicable to laboratory work. Please wear protective clothing, goggles and protective gloves when operating this apparatus.

Caution:

Instructed users only may operate the instrument.

Caution:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

Caution:

The On/Off switch does not disconnect the instrument from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.

Only set up the device on a solid, level surface.
Only use the device in a dry environment; do not use it in an explosive atmosphere.
If you want to move the device to another location, never hold it at the water level tube at the side, but at the condenser body and the lower part only.

Use only accessories and spare parts from the manufacturer. Use only the hoses supplied by the manufacturer. Do not use old hoses.

Repairs within the lower part of the device (electrical circuits, heating) may only be carried out by authorized electricians. Never send the device to your supplier for repair without first cleaning it.

4. Intended use

The distilling apparatus is used for mono-distillation of tap-water. It is designed to prepare approx. 4 litres of distilled water per hour.

5. Setup and operation

The apparatus consists of two parts: the glass part and the electrical part. Both parts are connected by a plastic flange. The upper part of the glass vessel contains the cooling coil (12). This upper part is separated by glass from the lower part. The vapour is passed through a tube, which is connected with the separator. By means of the cooling in the upper part, the vapour turns into distilled water, which runs off through the distillate drain (1) which is positioned above the glass separator.

5.1 Setup of the instrument

- Only set up the device on a solid, level surface.
- Please note that temperatures of up to 50 °C can be reached under the device.
- Only set up the device in a dry environment.
- The lower edge of the appliance must be horizontal and above the upper edge of the sink.
- The minimum safety distance from the water connection / water outlet should be 1 m to avoid the risk of splashing. Never place the device closer to it!

Caution:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

5.2 Installing a reservoir bottle

In order to achieve the conductivity of the distilled water as mentioned in the technical data, the apparatus and the reservoir bottle must form a closed system. Carbon dioxide from the air dissolves in water and would cause an increase in conductivity.

- Make sure that the reservoir bottle stands beneath the distillate drain (1).
- Place the reservoir bottle to the right of the apparatus.
- To connect the distillate drain (1) with the reservoir, use the silicone tubing.

5.3 Protective cap

Unscrew the protective cap (2) and place it on the lower part of the apparatus behind the level tube (3) for storage.

5.4 Connection of the tubing

Water supply: minimal pressure 200 kPa
 maximal pressure 600 kPa

Fluctuation of pressure: 10 % of water pressure max. In case of a higher rate of fluctuations, we recommend using a commercially available pressure regulator.

Use the PVC tubing supplied for the connections.

Before use immerse the end of the PVC tubing into hot water (of approx. 80 °C) for about 3 minutes. This will make the tubing flexible which now can be easily fitted on the glass tubing connections.

5.4.1 Connection of the cooling water

- Fasten the end of the tubing supplied (dia. ext. 19 mm) to the bent water connection (4) of the still body.
- Lead the tubing to the water tap and, if necessary, cut it off to suitable length.
- Fasten this end of the tubing to the water tap.

5.4.2 Connection of the water drainage

Water from the cooling coil (12) is led to the level tube (3) which also supplies the water to the evaporator. Surplus water is led over the overflow connection (5) to the sink. The water level in the evaporator is displayed at the point of the overflow connection (5) in the level tube (3).

- Fasten the end of the supplied PVC hose (dia. ext. 19 mm) to the overflow hose connection (5)
- Lead the other end of the tubing to the sink and cut it off to suitable length if necessary.
- Put this end into the sink and secure it against falling out.
- Do not allow the hose to sag.

Attention: When fitting the tubing to the overflow connection (5), support the level tube (3) with your other hand to avoid glass breakage.

5.4.3 Connection of the outlet valve (6)

Fasten one end of PVC tubing (dia. ext. 19 mm) to the outlet valve (6). Cut the other end off to suitable length and put it into the sink. Secure it against falling out!

5.5 Preparation of operation

Make sure that

- the power switch (7) is switched off (the green pilot lamp does not light up)
- the protective cap (2) is unscrewed
- the water tap is closed
- the stopcock (if any) of your reservoir bottle and the outlet valve (6) of the apparatus are closed

Now insert the power connector into the power socket.

5.6 Adjusting the water supply

- Turn on the water tap slowly to allow a small amount of water only to run through the overflow connection (5). The flow capacity then is approx. 0.5 to 1 litre per minute.
- As soon as the water in the evaporator (8) almost covers the heating coil, switch on the power switch (7). The green pilot lamp (9) lights up.
- As soon as the water in the evaporator (8) starts boiling, slowly reduce the water amount until vapour leaks out of the degassing outlet (10). Do not turn off the tap!
- Now slowly increase the water supply until it is just sufficient to stop the vapour from escaping. Now the water supply is properly adjusted.

Note: When operating the apparatus for the first time, please discard the water which has been distilled in the first hour of operation.

Attention:

Do not touch the still body! Danger of burns!

Attention:

Do not accidentally hit the glass vessel with hard objects. There is a risk of breakage and an additional risk of scalding while the device is in operation or shortly after it has been switched off.

Attention:

Risk of electric shock.

5.7 Interruption of water supply

If there is insufficient supplying of water or if the water flow is interrupted, vapour escapes from the degassing outlet (10). The water level in the evaporator declines. After approx. 5 minutes the heater will be switched off automatically.

Attention:

- Turn off the water tap
- Switch off the apparatus (7) and allow the apparatus to cool down for at least 15 minutes
- Correct the malfunction in the water supply.
- Switch off the temperature limit: take a piece of non-conducting material (length at least 50 mm and max. outside diameter 3 mm), insert it into the hole (11) on the reverse side of the lower part of the apparatus and press slightly inwards. Then remove the piece.
- To continue distilling, now repeat the steps described from 5.5.

5.8 To switch off

- Press the main switch (7). The green pilot lamp goes out.
- Close the water tap.
- Pull out mains supply plug out of the mains socket.

Attention: The water still must have cooled down completely before the remained water can be drained through the outlet valve (6). After that close the outlet valve.

5.9 Cleaning

Deposits (e.g. lime) can arise during the distillation process on the interior parts of the evaporator (8). This gradually reduces the quality of the distillate.

In case of considerable deposits in the evaporator or when the quality of the distillate no longer meets the requirements for the intended application, the evaporator must be cleaned. In case of heavy deposits in the cooling coil, they must be cleaned at a later time, however, these deposits do not affect the quality of the distillate.

The producer cannot accept any warranty claims if the heater is damaged as a result of overheating because the apparatus has not been sufficiently cleaned.

The following items are required for cleaning:

- protective clothing, gloves and goggles
- a suitable funnel
- a funnel stand, minimal height 70 cm
- a graduated jug of at least 2 litres capacity
- a suitable vessel for the disposal of the acetic acid
- acetic acid, 10 to 50 %
- tape water for rinsing

Warning:

Disconnect the apparatus from the mains by removing the plug from the socket before cleaning. There is a risk of corrosive burns. Wear protective clothing, gloves and goggles during the cleaning process. Dispose of the acetic acid in accordance with approved procedures. For your personal safety we recommend disconnecting the mains supply before cleaning. Keep your reservoir bottle connected. Allow the water still to cool down completely. We recommend to mark the apparatus with a label: "Do not switch on! Cleaning with acid."


5.9.1 Cleaning the evaporator

- Drain the water through the outlet valve (6) after the apparatus has cooled down completely.
- Close the valve (6)
- Now hang the other end of the tubing connected with the outlet valve (6) into the vessel for the disposal of the acetic acid. For safety reasons place the vessel into a collection tray.
- Firmly screw on the protective cap (2) on the pressure equalization opening of the level tube (3).
- Insert the funnel into the free end of the tubing leading to the overflow connection (5). If necessary, immerse the tubing into hot water (approx. 80 °C) for about 3 minutes before fitting. **Make sure the funnel fits tightly!**
- Now insert the funnel in the stand. The lower end of the funnel cone must be on the same level as the distillate drain (1) – see picture on the next page.
- Use the graduated jug to pour approx. one litre of acetic acid (10 to 50 %) into the funnel. Do not fill up the funnel with more than a third of its capacity to avoid spilling. Do not fill the evaporator to a higher level than the lower end of the distillate drain (1). You may regulate the filling speed as follows: To speed up the filling of the evaporator, place the funnel into a higher position. Place the funnel in a lower position to control the rate of addition of acid as it approaches the lower edge of the distillate drain tube (1), which is the point that must not be exceeded.

Attention:

Make sure that acid doesn't drip into the reservoir bottle. In case that acid already dripped into the reservoir bottle, carefully rinse the bottle and the bottle top.

- As soon as the deposits have dissolved, drain the acid from the evaporator via the outlet valve (6) into the vessel
- Close the outlet valve (6)
- Now hang the tubing of the outlet valve back into the sink.

Rinsing the evaporator

- Slowly pour approx. one litre of tap water into the funnel.
- Carefully fill the evaporator up to the lower edge of the distillate drain (1)
- Now drain the water through the outlet valve (6) and close the outlet valve again.
- Repeat the rinsing three times.
- Remove the funnel and hang the overflow tubing into the sink.


Unscrew the protective cap (2) and keep it on the base of the apparatus behind the level tube (3).

After the cleaning:

Let the apparatus run in distilling mode for one hour. Discard the distillate.

5.9.2 Cleaning of the cooling coil (12)

- After the apparatus has cooled down completely, drain the water out of the evaporator through the outlet valve (6). The outlet valve (6) must remain opened!
- Hang the tubing ends of the outlet valve (6) and the overflow connection (5) into the vessel for the disposal of the acid. For safety reasons place the vessel into a collection tray.
- Firmly screw on the protective cap (2) on the pressure equalization opening of the level tube (3).
- Disconnect the water supply tubing from the tap and insert the funnel into the free end of this tubing. **Make sure that the funnel fits tightly!**
- Now insert the funnel in the stand. The lower end of the funnel cone must be on the same level as the degassing outlet (1).
- Use the graduated jug to pour approx. 300 ml of acetic acid (10 to 50 %) into the funnel. Be carefully and avoid spilling!

Attention:

Acetic acid flows through the outlet valve and the overflow connection into the vessel.

- As soon as the deposits have dissolved, pour 300 ml of tap water into the funnel.
- After the liquid has flown through, hang the tubing of the outlet valve and the tubing of the overflow connection back into the sink.

Rinsing of the cooling coil

- Slowly pour approx. one litre of tap water into the funnel.
- As soon as the water has flowed through, close the outlet valve (6)
- Remove the funnel and reconnect the water supplying tubing with the water tap.
- Unscrew the protective cap (2) and keep it on the base of the apparatus behind the level tube (3).

After cleaning

Let the apparatus run in distilling mode for one hour. Discard the distillate.

6. Dismantling and transport

In case that the instrument has to be returned to the supplier for repair, please clean it carefully before dismantling.

6.1 Disassembly

- Switch off the instrument by means of the rocker switch.
- Disconnect the instrument from the mains.
- Now you can remove the instrument from the working area.
- Let the water run out of the apparatus completely.
- Use the supplied Allen key 3 mm to loosen all 8 screws of the clamping ring. Please open the screws step by step so that the pressure of the two-part clamping ring on the glass parts is released evenly. This will prevent breakage due to uneven pressure on the glass part.
- Loosen the screws completely.
- Lift the glass part with one hand and the clamping ring with the other.
- Then remove the silicone O-ring from the flange of the glass part.

6.2 Transport

Place the instrument and its parts in its original packaging or in a suitable container to protect it during transport. Close the packaging with adhesive tape.

In case you do not use the original packaging, please mark the box with the following notes:

- Glass symbol (Fragile! Handle with care!)
- Umbrella (Keep dry!)
- Content (List of content)

6.3 Reassembly

- Insert one silicone O-ring into the groove of the base.
- Then insert the second silicone O-ring into the flange of the glass part.
- Now place the flange of the glass part in the center of the base. The flange must completely fit into the recess of the base.
- Now place the two halves of the clamping ring on top of the base. The holes for the screws on the clamping ring must be exactly above the holes in the base. The second silicone O-ring must fit into the groove of the clamping ring.
- Now insert the screws and use the supplied Allen key 3 mm to tighten the screws step by step. Never fully tighten one half of the clamping ring while the other half remains unscrewed. This will prevent breakage due to uneven pressure on the glass part.
- Tighten the screws until the distance between the lower edge of the clamping ring and the base part is around 1 mm.

7. Disposal

Damaged and/or waste electric or electronic equipment have to be disposed of at the designated Recycling Depot. Also the packing material should be disposed of environment-friendly (material separation)

8. Warranty and Liability

Paul Marienfeld GmbH & Co. KG provides a guarantee for a period of 24 months from date of delivery for the instrument purchased by you. This warranty applies only to defects in material or workmanship. There is no warranty for defects or malfunctions caused by neglect, improper use or improper service.

In case of defects in material or workmanship, the instrument will be repaired or replaced. In case of replacement, the warranty for the new unit will persist for the remaining period of the primary warranty period. Charges for returning the instrument have to be borne by the customer.

Any claims regarding guarantee and liability are limited to professional and proper use of the instrument and observance of the general safety instructions and this instruction manual.

Paul Marienfeld GmbH & Co. KG is not liable for any material or personal damages caused by improper use and does not accept liability for consequential damages. The warranty expires in case of removal or alteration of the identification label and/or the serial number of the instrument.

9. Technical Data

Mains:	230 V, 50 Hz
Power consumption:	approx. 3100 W
Nominal current:	approx. 14 A
Electrical protection – tube fuse 5 x 20 F 8 A:	1 A
Output of distillate:	approx. 4 l/h
Required consumption of cooling water:	approx. 40 l/h
Temperature of distillate:	approx. 95 °C
Conductivity of distillate at 25 °C:	1,5-2,5 µS/cm
Dimensions of apparatus:	190 x 190 x 600 mm

10. Declaration of conformity

We hereby declare that this product

Distilling apparatus for water DP 4000

is in conformity with the following directives:

EMV directive 89/336/EEC amended by 91/263/EEC
EMV directive 72/23/EEC

In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.

Paul Marienfeld GmbH & Co. KG
Harry Marienfeld, President¹

¹ The above declaration has been generated automatically and is valid without signature

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